

Molecular Collection Device/Test Guide

Cytology Primary Collection Devices



ThinPrep Pap

Test (Stability):

Pap (6 weeks)
 HPV (14 days)
 HPV Genotyping (14 Days)
 CT/NG (30 Days)
 Trichomonas (14 Days)
 Dual Stain: p16/Ki67



Sterile Specimen Cup

Test (Stability):

Urine for Cytology (48 Hours)
 CT/NG (24 Hours)
 M gen (24 Hours)
 Trichomonas (24 Hours)

Molecular Primary Collection Devices



Aptima® Multitest Swab Kit

Test (Stability):

CT/NG (60 Days)
 Trichomonas (60 Days)
 Candida Vaginitis/Trichomonas (30 Days)
 Bacterial Vaginosis (30 Days)
 HSV 1&2 (36 Days)
 M gen (60 Days)



Aptima® Urine Collection Kit

Test (Stability):

CT/NG (30 Days)
 M gen (60 Days)



eSwab (Liquid Amies)

Test(Stability):

GBS (24 Hour)